

WHITE PAPER

SEPT 2025

The Mock-Up: An Ultimate Guide to Mastering All 3D Printing Options

This white paper explores the evolution of mock-up fabrication in digital dentistry, detailing three 3D printing protocols that cater to different clinical needs: a printed model used to produce a PVS index, a direct printed index with a biocompatible resin, and a printed mock-up with tooth-colored resins. This guide outlines the full workflow for each and provides dental professionals with the guidance to choose the most suitable method for their specific cases.

Contents

About This Guide

4

The Concept Around the Mock-up

5

Types of Mock-ups

5

Motivational Mock-up: Creating the Emotional Connection With the Patient

5

Technical Mock-up: Technical Decisions and Clinical Execution

6

The Evolution of the Mock-up: From Past To Present

7

Direct Mock-up Technique

7

Indirect/Direct Mock-up Technique

8

Fully Indirect Mock-up Technique

9

Modern Mock-up Options Using 3D Printing

10

3D Printing Empowered Mock-up Techniques

10

Overview of 3D Printing Protocols for the Production of the Mock-up

10

Decision Factors for Choosing the 3D Printing Approach

11

Advantages and Challenges of Each Protocol

12

Protocol 1: Model Production Used To Fabricate a PVS Index

13

1. Design Considerations

14

2. Material Selection

14

3. Nesting in PreForm Dental

17

What Would Be the Ideal Setup?

20

4. 3D Printing and Post-Processing

21

5. Clinical Execution

22

Protocol 2: Printed Index With IBT Flex Resin

24

1. Design Considerations

25

2. Material Selection

26

3. Nesting

27

4. 3D Printing and Post-Processing

27

5. Clinical Execution

29

How does this material and protocol compare to a traditional PVS produced index?

30

Protocol 3: Printed Mock-up Shell With Tooth-Colored Resins**31**

1. Design Considerations 32

2. Material Selection 32

3. Nesting in PreForm Dental 34

4. 3D Printing and Post-Processing 35

5. Clinical Execution 37

Conclusion 39

About This Guide

[Formlabs Dental](#) and [Digital Smile Design \(DSD\)](#) have joined forces to redefine what's possible in modern dentistry. Since 2019, this collaboration has brought together DSD's pioneering vision in comprehensive guided dentistry and Formlabs Dental's leadership in [advanced dental 3D printing technologies](#). By combining clinical innovation with cutting-edge manufacturing tools, we've been working to develop the next generation of digital workflows and practical solutions for dental professionals.

As part of this joint effort, we've explored and tested a variety of protocols to create the best-in-class mock-ups, essential tools for diagnostic, clinical, and communication purposes. DSD has continuously pushed the boundaries of what 3D printing can offer in the dental space, providing invaluable feedback that has shaped product development and inspired new clinical applications. Meanwhile, Formlabs Dental has delivered the [materials](#), hardware, and [reliability](#) to turn these ideas into tangible, everyday solutions.

Together, we have witnessed the significant evolution of the mock-up technique. What once began with printing physical models to create PVS indexes has now advanced to the direct printing of biocompatible, end-use appliances. Today, 3D printing offers incredibly flexible approaches, streamlining workflows and enhancing patient outcomes.

This guide presents three proven protocols for using 3D printing to fabricate mock-ups, offering flexible, efficient, and innovative options to suit a wide range of clinical needs. Below, we explain:

- How to choose the best 3D printing workflow depending on the type of mock-up you need to produce
- Provide guidance on selecting the right approach
- Walk you through the full 3D printing process from design to final outcome


This comprehensive guide is the result of many years of collaborative work. It has been prepared after multiple intensive sessions, extensive testing, and valuable insights gained from practical application, including during DSD Clinical Courses. The authors, Dr. Luken de Arbeloa and Dr. Elisa Praderi, have dedicated years to bringing these concepts together into this consolidated resource.

The Concept Around the Mock-up

The mock-up has long been a cornerstone of restorative and esthetic dentistry, a trusted tool to help patients visualize potential outcomes. In recent years, DSD has redefined its role. By integrating psychology, 3D design, and structured clinical workflows, DSD transformed the mock-up into a powerful bridge between patient motivation and clinical execution. No longer just a visual aid, the mock-up is now a communication tool, a motivational driver, and a technical reference — all depending on the timing and intent within the treatment journey.

At the core of the DSD approach is the understanding that mock-ups serve two very different purposes: **motivation and technical execution**. These are separated by a key clinical moment: **case acceptance**.

If we look at the DSD Workflow, the listed steps are followed:



Based on the above, we can identify that we would use different types of mock-ups for different phases and purposes.


Types of Mock-ups

Motivational Mock-up: Creating the Emotional Connection With the Patient

The motivational mock-up is designed to inspire and engage the patient before case acceptance. It serves as an emotional tool to help the patient visualize their future smile, build trust, and understand the value of the proposed treatment. At this stage, clinical accuracy is not the priority as we are strictly focusing on motivation only.



From a design and 3D printing perspective, the designed mock-up is intentionally bulky and fully additive, with no preparation guides or fine anatomical detail, as it isn't intended to support technical decisions. The typical workflow involves printing a model used to fabricate a PVS index and using bis-acryl resin to transfer that design to the patient's mouth.




Since this mock-up is used primarily for photos or social-distance visualization, ultra-high precision isn't required — layer heights of **100 to 200 µm** are more than adequate. This makes it ideal for **fast prints using cost-effective materials**. For example, we can print a model in 160 or 200 µm layer height with [Form 4B](#) in less than 10 mins.


Another alternative that has appeared most recently is the use of biocompatible tooth-colored resins for the production of 'the shell'. For example, we can print the motivational mock-up using [Premium Teeth Resin](#) at 100 µm that will take 20 minutes. This process becomes more efficient as there are less steps involved and less material consumption.

Technical Mock-up: Technical Decisions and Clinical Execution

Once the patient has accepted the treatment plan, the focus shifts to the technical mock-up, which plays a crucial role in **supporting clinical decision-making and ensuring precise treatment execution**.



Unlike the motivational mock-up, the technical version is more refined and functional, featuring accurate anatomical contours that closely match the planned final result. It guides important clinical steps such as crown lengthening, tooth preparation depth, and soft tissue management.



The workflow may include **direct printed shells** or **detailed mock-ups transferred using printed indexes**. Precision is critical at this stage, requiring high-resolution 3D printing and careful material selection to ensure reliable handling and fit. Ultimately, the technical mock-up validates and guides the clinical protocol, helping dental professionals execute complex restorative procedures with confidence, consistency and predictability.

The following table summarizes the concepts of the different types of mock-ups:

	MOTIVATIONAL MOCK-UP	TECHNICAL MOCK-UP
The moment	Before case acceptance	After case acceptance
The design	Adapted design to allow placement in the mouth over the existing pre-op situation	Ideal design to guide clinical procedures and the final outcome
The purpose	- Emotional purpose - Motivate the patient to do the treatment - Wow effect	- Rational and technical purpose - Technical communication with patients and/or team - Guide clinical procedures
The material	Bis-acryl or shell	Bis-acryl or shell

The Evolution of the Mock-up: From Past To Present

The DSD team has already prepared valuable guides that offer a thorough analysis of the mock-up's concept and the evolution of the technique. For more detailed information, you can refer to "[The Complete Guide to the DSD Mockup Success Formula](#)."

Historically Used Protocols

DIRECT MOCK-UP TECHNIQUE		INDIRECT/DIRECT MOCK-UP TECHNIQUE		FULLY INDIRECT MOCK-UP TECHNIQUE		
Fully direct manual approach	Prefabricated direct composite shells	Indirect/direct manual	DSD Mock-up	Fully indirect handmade shell	Fully indirect digital snap-on	DSD Shell

Direct Mock-up Technique

1. Fully Direct Manual Approach: where composite is used directly in the patient's mouth to do the smile design. The challenge here is that the clinician must have a great sense of facial esthetics. This step can be time-consuming and labor-intensive, and it may be challenging to convey its value to the patient and justify the time investment at this stage.



Case by Dr. Galip Gurel
From: "The Complete Guide to the DSD Mockup Success Formula"


2. Prefabricated Direct Composite Shells: Historically, one of the key analog approaches was the direct mock-up technique, developed by Dr. Paulo Kano. This method uses prefabricated composite shells that are positioned directly in the patient's mouth. The process involves manual waxing upon a model, followed by copying the model's texture using silicone PVS index to create these composite shells. While effective, this technique is very time-consuming and requires considerable manual skill to execute properly.



Case by Dr. Paulo Kano
From: "The Complete Guide to the DSD Mockup Success Formula"

Indirect/Direct Mock-up Technique

1. Indirect/Direct Manual: An initial impression of the current patient situation is needed, followed by a manual wax-up on a model that may be the product of the measurements from a 2D DSD Design. The waxed-up model must be transferred to the patient's mouth using a PVS silicone index filled with bis-acryl resin. The shape and texture in these cases is applied manually by the dental technician imitating the anatomy of the teeth and the aesthetic result depends on the correct interpretation of the 2D DSD Design.



Case by Dr. Christian Coachman

From: "The Complete Guide to the DSD Mockup Success Formula"

2. DSD Mock-up: This is an evolution of the indirect/direct manual technique, substituting the manual wax-up for a digital one. The main benefits are:

- The creation of the Motivational Mock-up or wax-up after completing a DSD Ideal Design
- The overlapping of the .STL files with the patient's face and real tooth textures

As there is an intermediate digital step, the mock-up is done in the second appointment.

To create the mock-up, the design is 3D printed as a model and a silicone index or vacuum tray created from it. The way to bring this design to the patient's mouth is exactly the same as the previous mock-up technique using a PVS silicone index, which in this case is a copy of the Motivational 3D Printed Model using bis-acryl.



Case by Dr. Ricardo Brito

From: "The Complete Guide to the DSD Mockup Success Formula"

Fully Indirect Mock-up Technique

1. Fully Indirect Handmade Shell: This mock-up is fully fabricated in the dental lab based on a handmade wax-up and without the involvement of a patient. From this, removable acrylic or bis-acrylic shells are produced manually on a model by copying the wax-up. The device requires adequate thickness and connector size to resist the path of insertion without breaking.



Wax-up and mock-up performed by Well Lab, São Paulo

From: "The Complete Guide to the DSD Mockup Success Formula"

2. Fully Indirect Digital Snap-On: This mock-up is achieved using a digital wax-up. The project is adapted in CAD/CAM software giving it the same characteristics as the fully indirect handmade shell including adequate thickness and strong connectors to resist the path of insertion. This usually makes these devices more bulky, less delicate, and less aesthetic. However, it can be used as an orthopaedic/functional device to equilibrate the bite.



Wax-up and mock-up performed by Well Lab, São Paulo

From: "The Complete Guide to the DSD Mockup Success Formula"

3. DSD Shell: Due to the evolution of technology, this fully indirect mock-up technique is different from the traditional DSD Mock-up technique in the way the design is transferred to the patient's mouth.

Taking advantage of 3D printing technology, a biocompatible tooth-colored resin is used to print the mock-up. The design is created by modifying the 3D DSD Ideal Design into a Motivational Design. This 3D printed mock-up is ready to be placed in the patient's mouth without the need for a silicone index. It consists of three very fine sections that are then adhered to the teeth with flowable composite, obtaining a result that is facially integrated with natural textures and adequate color, and which is quick and easy to place.



DSD Shell mock-up by Dr. Luken de Arbeloa

From: "The Complete Guide to the DSD Mockup Success Formula"

Modern Mock-up Options Using 3D Printing

3D Printing Empowered Mock-up Techniques

INDIRECT/DIRECT MOCK-UP TECHNIQUE

THE DSD MOCK-UP



In the current landscape, the DSD mock-up represents a hybrid analog-digital method. After designing the smile digitally, the model is 3D printed, and a PVS silicone index is created for transferring the mock-up intraorally with bis-acryl resin. This technique combines the precision of digital design with the familiar analog transfer process.


PRINTED INDEX



A more recent innovation is printing the transfer index directly using a flexible biocompatible material (eg. [IBT Flex Resin](#)). This eliminates the need to print a physical model, streamlining the workflow by saving both time and materials. This method retains the clinical effectiveness of the traditional transfer while introducing digital efficiency.

FULLY INDIRECT MOCK-UP

DSD SHELLS



The fully indirect mock-up has embraced digital fabrication through the production of 3D printed shells using biocompatible, tooth-colored resins. These digitally fabricated mock-ups offer excellent aesthetics, precision, and cost efficiency. Their quality and convenience have made them a favorite among clinicians looking for a completely digital solution to mock-up fabrication.

Overview of 3D Printing Protocols for the Production of the Mock-up



Protocol 1: Model Production Used to Fabricate a PVS Index

The foundational DSD mock-up. A printed model which is used to create a PVS index, transferred to the mouth with bis-acryl resin.



Protocol 2: Printed Index With IBT Flex Resin

Skip the model entirely. Print a flexible index directly, simplifying the workflow and reducing material use.




Protocol 3: Printed Mock-up Shell With Biocompatible Tooth-Colored Resins

A fully digital, totally indirect mock-up. High-precision, aesthetic, and ideal for chairside try-ins or as a mid-term provisional.

From a 3D printing perspective, each method has a specific recommended material and approach. Motivational mock-ups require less precision (layer heights of 100–200 μm are sufficient), while technical mock-ups demand maximum accuracy and material fidelity (50 μm).

By understanding the purpose behind each type of mock-up and applying the right 3D printing strategy, dental professionals can elevate both patient engagement and clinical outcomes, combining emotion and precision like never before.

In the **full digital workflow utilizing 3D printing**, the following steps are covered:



Decision Factors for Choosing the 3D Printing Approach

In addition to connecting the motivational and technical mock-ups to the right 3D printing materials and strategies, there are other questions that dental professionals need to ask themselves to confirm they are in the right position to adopt each workflow.

	PROTOCOL 1: MODEL PRODUCTION USED TO FABRICATE A PVS INDEX	PROTOCOL 2: 3D PRINTED INDEX WITH IBT FLEX RESIN	PROTOCOL 3: PRINTED MOCK-UP SHELL WITH TOOTH-COLORED RESIN
Material classification	Non-biocompatible resin	Biocompatible resin Class I (US & EU)	Biocompatible resin Class II (US) & Class IIa (EU)
Type of mock-up	Motivational	Technical	Motivational Technical
Resin regional availability	Available globally	Available in specific regions	Available in specific regions
Production time	Short: print a model in less than 10 min	Long: takes the longest to print	Intermediate: can print in ~20 min
Efficiency level	Intermediate: requires time to print the model and fabricate the traditional silicone index upon the model	Low: requires time to design the index in the CAD software	High: printed shell can be used intraorally
CAD steps involved	1. Digital wax-up	1. Digital wax-up 2. Index design upon wax-up	1. Digital wax-up 2. Shell design
Level of CAD expertise	Intermediate	High	Intermediate

	PROTOCOL 1: MODEL PRODUCTION USED TO FABRICATE A PVS INDEX	PROTOCOL 2: 3D PRINTED INDEX WITH IBT FLEX RESIN	PROTOCOL 3: PRINTED MOCK-UP SHELL WITH TOOTH-COLORED RESIN
Protocol level of complexity	Low	High	Intermediate
Material observations	N/A	• When printing with biocompatible resins, it's mandatory to follow the Instructions for Use and Manufacturing Guide for each material. • Biocompatible resins require the use of a dedicated Build Platform and Form Wash. • Following the exact wash times, dry times, and post-curing times and temperatures is fundamental to guaranteeing success of the printed part. • Biocompatible resins come with specific clinical indications and minimum thickness requirements. Deviating from these guidelines can lead to off-label use of the materials, potentially exposing patients to increased risks.	

Advantages and Challenges of Each Protocol

	PROTOCOL 1: MODEL PRODUCTION USED TO FABRICATE A PVS INDEX	PROTOCOL 2: 3D PRINTED INDEX WITH IBT FLEX RESIN	PROTOCOL 3: PRINTED MOCK-UP SHELL WITH TOOTH-COLORED RESIN
Advantages	• Material availability is not a problem. • Fastest option for motivational mock-up. • Requires one CAD step. • Best protocol for beginner users.	• Interproximal reproducibility and transfer of natural textures • Great option for a motivational mock-up and a DSD Direct technique. • Best protocol for advanced users. • A separating agent might not be needed when using the printed index with bis-acryl resin- users shall confirm material compatibility extraorally.	• Directly print the mock-up. • The most efficient workflow with minimum material consumption. • Best protocol for beginner to advanced users.
Challenges	• Requires a combination of analog and digital steps. • Producing two parts means more material consumption. • The PVS index, though a proven method, requires care while creating it and performing the cuts.	• Regional availability of material. • Requires two CAD steps. • CAD software limitation. • If using a traditional or flowable composite, it might require the use of a separating agent. • Technique-sensitive: post-processing protocol is key.	• Regional availability of material. • Requires two CAD steps. • Material has a minimum thickness recommendation. • Technique-sensitive: post-processing protocol is key.

Protocol 1: Model Production Used To Fabricate a PVS Index



The first option is to print the wax-up model, which professionals can then use to fabricate the PVS index.

Tools needed to do this protocol:

DIGITAL CAPTURING AND DESIGN	DIGITAL DESIGN	DIGITAL PRODUCTION	CLINICAL TOOLS
Intraoral scanner Digital camera	CAD software for doing the smile design	Form 4B 3D printer Resin Tank Resin Mixer Flex Build Platform Form Wash ^(2nd Generation) Form Cure ^(2nd Generation) Fast Model Resin Precision Model Resin Accessories: IPA (99%) Protective equipment	PVS Scalpel Bis-acryl resin Probe Tweezers Gloves Black marker Pink composite (if necessary)

Workflow Steps

FIRST APPOINTMENT: DATA ACQUISITION

Intraoral scanning
Photography
Videos

PLAN AND DESIGN

Design the mock-up and export the model with a wax-up for printing

3D PRINT AND POST-PROCESS THE PRINTED MODEL

Material selection
Nesting
3D Printing
Post-processing

FABRICATE THE INDEX

Use the model to fabricate the PVS index and cut the index

CLINICAL EXECUTION

Use PVS index to do the mock-up using bis-acryl resin

1. Design Considerations

- The wax-up will be done on the imported data that was acquired in the first appointment with the intraoral scanner and digital pictures and videos.
- For the **motivational mock-up**, the wax-up design should be only additive and designed mainly in the buccal area.
- For the **technical mock-up**, the wax-up is done following the traditional principles, looking for proper function and performance.
- There are different softwares available that you can use to do the digital wax-ups such as [3shape](#), [exocad](#), [NemoTec](#), among others.

2. Material Selection

When analyzing the 3D printing workflow, one of the important steps is selecting which model resin to use.

There are different model materials available from Formlabs Dental which can fulfill this application including [Fast Model Resin](#), [Grey Resin V5](#), and [Precision Model Resin](#). In order to select the right material, it's important to analyze the requirements for each application.

- For **motivational mock-ups**, we need a **fast solution** as the objective is printing a model with the wax-up that doesn't need to be the most precise option. Based on this criteria, the best option is to use Fast Model Resin. This material enables the printing of models in 100, 160, and 200 µm layer height settings. Any of these layer height options work for motivational mock-ups. It's up to the dental professional to define the layer height based on their personal preferences.
- For **technical mock-ups**, we need a solution that provides a **smooth surface finish and precision** so we can accurately transfer the planned design to the patient's mouth. In this case, Precision Model Resin printed at a 50 µm layer height setting would offer the best surface accuracy. The PVS indexes produced upon these models will accurately transfer the anatomical features and textures designed.

Though there are other materials available such as Grey Resin V5 and [White Resin V5](#), these are less versatile for other dental applications. Additionally, neither offers the speed of Fast Model Resin nor the precision of Precision Model Resin.

The following table summarizes the main differentiating points between Fast Model Resin and Precision Model Resin.



FORMLABS RESIN	FAST MODEL RESIN	PRECISION MODEL RESIN
Ideal for	Printing a model to create a PVS index for a motivational mock-up.	Printing a model to create a PVS index for a technical mock-up.
Other uses of the material	• Diagnostic models • Models for thermoforming clear aligners, retainers or bleaching trays • Models used for fitment check of occlusal splints or hawley retainers • Antagonist models	• Removable die models • Implant analog models • Restorative models with soft tissue • Models with wax-up to fabricate clear silicone indexes for injectable or pressed composite technique
Color	Light Grey	Beige
Cost of cartridge	\$/€ 79	\$/€ 129
Cost per part	\$/€ 1-3	\$/€ 2-4
Layer height recommended and estimated print time for one model	Wax-up model for motivational mock-up: 200 µm: 8 min* 160 µm: 9 min* Other layer height option: 100 µm: 22 min* <i>*Please note that print times may vary depending on the geometry and height of the models</i>	Wax-up model for technical mock-up: 50 µm: 49 min* Other layer height options: 100 µm: 20 min* 25 µm: 2 h 45 min*
Post-processing time	Wash: • Form Wash (2nd Generation) at High Speed using IPA (99%): 5 min Dry Cure (choose one): • Form Cure (2nd Generation): 1 min, no heat • Form Cure (1st Generation)/ Form Cure L: 5 min, no heat	Wash: • Form Wash (2nd Generation) at High Speed using IPA (99%): 5 min Dry Cure (choose one): • Form Cure (2nd Generation): 5 min, no heat • Form Cure (1st Generation)/ Form Cure L: 5 min, 35 °C



FORMLABS RESIN

FAST MODEL RESIN

PRECISION MODEL RESIN

What type of mock-up can I use this for?

Motivational mock-up:
You can print a model in as little as 8 minutes. Use any of the layer heights for this purpose. For fastest results, use 200 µm or 160 µm.

Motivational mock-up:
Use 100 µm
Technical mock-up:
Use 50 µm

Guidance for choosing the right resin

If speed is important, this should be your go-to resin.

If you are looking to get the highest precision and accuracy, this material is the best choice.

Don't choose this material if speed is fundamental for you and print time is a concern.

If you are looking to have one resin that can fulfill the needs of highly accurate restorative models as well as printing models for doing the wax-up, then you can choose this resin.

When considering the motivational mock-up, Fast Model Resin is the ideal choice as it combines the fastest print speed, 160 or 200 µm layer height setting, and reduced material cost.



TYPE	SOLID: MODEL IS FILLED	HOLLOW: MODEL IS EMPTY	HONEYCOMB: MODEL HAS A HONEYCOMB STRUCTURE
Pros	Safe to print directly on the build platform, which will also be the fastest printing option.	Requires less resin to print the part.	Requires less resin to print the part compared to solid models.
	The base height can be reduced as the model is completely filled.	Can be printed directly on the build platform if the model has a minimum wall thickness of 3 mm and drain holes.	
	No drain holes needed.		
Cons	Some extra resin is used to print the interior part of the model.	More prone to deformations if minimum wall thickness is not followed.	
		Higher risk of print failure if no drain holes are added when printing models directly to the build platform.	
Observations	This is the preferred design as the model is completely filled and uniform. It's also self-supported.	Requires a minimum wall thickness of 3 mm.	
		Drain holes are required if printing the model directly on the build platform.	


3.2 Print Orientation

Print orientation is one of the variables that affects print time, the number of models that can fit in one build, and the direction the layers are built. When printing models, there are two possible options:

- Printing the models flat against the build platform with their base directly on the build platform surface.
- Printing models vertically or angled with support structures.

Each option offers different advantages and disadvantages, which are summarized below.

		
PRINT ORIENTATION	FLAT AGAINST THE BUILD PLATFORM	VERTICAL OR ANGLED
Pros	Fastest to print (under 10 min) Self-supported model: No supports are needed, less resin consumed Very accurate	Takes longer to print (35 min) Requires support structures to print the models, more resin used for these structures
Cons	Limits the amount of models you can fit on one build platform (11 models)	Maximizes the space and allows combining several models (21 models)



A full build platform of horizontally oriented models (top) versus vertically oriented models (bottom).

3.3 Layer Height

Fast Model Resin and Precision Model Resin offer different layer height options, as illustrated in the comparative table. Clinicians can select the material based on the mock-up type, desired print time, or cost per part, while also considering the optimal layer height. It's important to remember that layer height directly influences print time and surface finish. The table below provides a summary of these concepts.


LAYER HEIGHT	50 μm	100 μm	160 μm	200 μm
Mock-up type	Technical	Motivational	Motivational	Motivational
Print speed	Slowest	Intermediate	Fast	Fastest
Surface finish	Smoothest surface finish Excellent for accurate transfer of all details	Smooth surface	Smooth surface with some visible layer lines	Less smooth surface with some visible layer lines but works for a motivational mock-up
Available with Precision Model Resin	✓	✓	–	–
Available with Fast Model Resin	–	✓	✓	✓

What Would Be the Ideal Setup?

For both **motivational mock-ups** and **technical mock-ups**, the best recommendation is to print models that are **solid** and **placed with their base directly on the build platform**. The layer height selection would be dependent on the type of mock-up needed and professional preference.

1. Ideal Set-Up for Motivational Mock-up:

- Parameters
- Fast Model Resin
 - 160 µm
 - Solid model
 - Printed directly on the Build Platform



Printer
 **AbleKiwi**
Form 4B

Status

Not ready

Connection

Local

Cartridge

Fast Model V1

Tank

Fast Model V1

Job Settings

Material

Fast Model V1

Layer Thickness

0.160 mm

Print Settings

Fast Arch

Summary ⓘ 

 Total Print Time

8 m

 Volume

8.66 mL


 Layers

75

Print time and volume

2. Ideal Set-Up for Technical Mock-up:

- Parameters
- Precision Model Resin
 - 50 µm
 - Solid model
 - Printed directly on the Build Platform



Printer
 **CrispDingo**
Form 4B

Status

Finished

Connection

Local

Cartridge

Precision Model V1

Tank

Precision Model V1

Job Settings

Material

Precision Model V1

Layer Thickness

0.050 mm

Print Settings

Default

Summary ⓘ 

 Total Print Time

46 m

 Volume

8.67 mL

 Layers

240

Print time and volume

4. 3D Printing and Post-Processing

Following material selection and nesting in [PreForm Dental](#), the parts are printed and post-processed. The 3D printing workflow involves several steps: printer setup, washing, drying, and post-curing. The total production time, including these final steps, varies per material used.

4.1 Setting up Form 4B

When printing with the Form 4B 3D printer, the following are required: a resin tank, resin mixer, build platform, and resin cartridge of the selected material. Once these are correctly inserted the print job can be sent from PreForm Dental to the printer.

4.2 Post-processing

After the print job is completed, models must be post-processed.

First, parts are removed from the Flex Build Platform using the quick release technology. This enables users to streamline part removal and avoid breaking or damaging models.

Parts are then washed in IPA (99%) at high speed in the Form Wash (2nd Generation) and then dried with compressed air. Lastly, parts are cured in the Form Cure for a specific time and temperature (see the material instructions). The following graphics summarize the overall process and final production time.



Motivational Mock-up: Fast Model Resin at 160 µm



Print directly on the Build Platform




Print time:

1 model
160 µm
8 min



Wash time:

5 min




Cure time:

Form Cure: 5 min, no heat
OR
Form Cure (2nd Gen): 1 min, no heat



15-18 min production time


Technical Mock-up: Precision Model Resin at 50 μ m




>1 h production time

5. Clinical Execution

The silicone (PVS) index can be created directly on the printed model. There is no need to apply any separating agent to the model. It's important to bear in mind the thicknesses of the PVS material in each section. Also, the cuts of the index need to be properly done as shown in the following sequence.



Once the index is created, it will be used to do the mock-up in the patient's mouth. The full clinical procedure can be seen in the "[The Complete Guide to the DSD Mockup Success Formula](#)".



Protocol 2: Printed Index With IBT Flex Resin




Another 3D printing option that works for both motivational and technical mock-ups is to print the index with a biocompatible resin to transfer the designed wax-up to do the mock-up in the patient’s mouth.

Tools needed to do this protocol:

DIGITAL CAPTURING AND DESIGN	DIGITAL DESIGN	DIGITAL PRODUCTION	CLINICAL TOOLS
Intraoral scanner Digital camera	CAD software for doing the smile design and to design the index	Form 4B 3D printer Resin Tank Resin Mixer Flex Build Platform Form Wash (2nd Generation) Form Cure (2nd Generation) IBT Flex Resin Accessories: IPA (99%) Protective equipment	Bis-acryl resin Probe Tweezers Gloves Black marker Pink composite (if necessary)




Workflow Steps



1. Design Considerations

In this case, follow similar steps to those in protocol one during the design phase. Do the wax-up on the intraoral data that was captured from the patient. After the wax-up is done, the design of the index can be done upon it. Important parameters to bear in mind while doing the design of the digital index:

- **Minimum thickness:** The walls of the index need to be 4-5 mm thick.
- **Limits:** Adjust the design to fit perfectly in the cervical area.
- **Undercuts and interproximal areas:** Observe the undercuts you have and need. This technique enables replicating undercuts in a very reliable way.
- **Offset from teeth:** This needs to be zero.
- **Flat plane:** The part of the index that will be facing toward occlusal needs to be designed with a flat surface so it can be printed directly to the build platform.
- **Extension:** For motivational mock-ups, design a full-arch index. Define the palatal cut in a precise way to facilitate the cleaning process (removal of excesses of the bis-acryl resin).



2. Material Selection


The material indicated for this purpose is [IBT Flex Resin](#), a biocompatible resin used to print flexible parts for short-term use. The main properties of the material can be found in the following table:



FORMLABS RESIN	IBT FLEX RESIN
Ideal for	Printing the index for motivational or technical mock-up.
Other uses of the material	A Class I biocompatible resin for producing appliances for short term use: • Direct composite restoration guides for injectable or pressed composite techniques • Indirect bonding trays for bracket placement
Color	Clear
Cost of cartridge	\$/€ 249/L
Cost per part	\$/€ 3-5
Layer height recommended and estimated print time for one full-arch index	Motivational mock-up: 100 µm: 1 h 7 min* Technical mock-up: 50 µm: 3 h 28 min* <i>*Please note that print times may vary depending on the geometry and height of the parts</i>
Post-processing time	Wash: • Form Wash (2nd Generation) at High Speed using IPA (99%): 20 min • Then either rinse parts completely with fresh IPA from a spray bottle, or soak parts in fresh IPA (99%) for 10 min Dry (allow parts to dry for 30 min before curing) Post-cure: Submerge parts in a transparent container filled with water and place inside the cure unit. Cure units that you can use and their parameters are: • Form Cure (2nd Generation): 7 min at 70 °C • Form Cure (1st Generation)/Form Cure L: 30 min at 70 °C
What type of mock-up can I use this for?	Motivational mock-up: 100 µm Provides the fastest print with this material but will have more visible layer lines than 50 µm Technical mock-up: 50 µm Provides the best surface finish for the most accurate transfer of the designed textures intraorally
Guidance for choosing the right resin	Dental professionals should verify regional availability of this material here .

3. Nesting

When designing one of the planes of the tray to be flat, the parts shall be placed with this surface directly on the build platform. This will provide the fastest print speed as the model is self-supported.



4. 3D Printing and Post-Processing

Following nesting in PreForm Dental, the parts are printed and post-processed. The 3D printing workflow involves several steps: printer setup, printing, washing, drying, and post-curing. Total production time varies depending on the layer height chosen for printing.

4.1 Setting up Form 4B

Insert a build platform, resin tank, resin mixer, and IBT Flex Resin cartridge in the printer. Once the printer is set up, send the print job from PreForm Dental to the printer.


When printing with biocompatible resins, a dedicated build platform, resin tank, and resin mixer are needed to ensure biocompatibility.

4.2 Post-processing

Post-processing of parts printed with biocompatible resins requires following the steps established in the [Manufacturing Guide](#).

Wash:

- After parts are removed from the build platform, they need to be washed in IPA (99%) in the Form Wash (2nd Generation) at High Speed in two steps (20 min + 10 min). It's important to use a dedicated Form Wash for washing parts printed with IBT Flex Resin (the wash should not be used for parts printed in other materials) and follow the exact timing as presented in the [Manufacturing Guide](#).
- Immediately after washing, it's important to dry the printed parts with compressed air to ensure that there is no remaining IPA or resin on the printed part, especially in the intaglio surface. This is critical to ensure proper fit of the appliance during the clinical execution. The following picture shows how an insufficient wash can alter the intaglio surface.
- Allow parts to dry for at least 30 minutes before proceeding to post-curing.




Cure: Parts need to be finally cured in the Form Cure (1st Generation) or Form Cure (2nd Generation) for the time and at the temperature in the IBT Flex Resin [Manufacturing Guide](#). This final step will ensure parts are ready to be used. The following graphics summarize the overall process and production time.

Motivational Mock-up: IBT Flex Resin at 100 μ m



Technical Mock-up: IBT Flex Resin at 50 μ m



5. Clinical Execution

In this case, the index is going to be used just like a traditional PVS index with bis-acryl resin. No separating agent is required when using bis-acryl resin with the printed index. However, users should always confirm extraorally that their chosen bis-acryl material does not bond to the printed parts before intraoral application.




How does this material and protocol compare to a traditional PVS index?

In the digital workflow, we have full control over the thickness of the material, allowing us to create precise, efficient designs in just 15–20 minutes. This not only improves accuracy but also leads to better clinical results with minimal excess. One key advantage is the precision in the palatal area, making it easier to clean and manage during placement. Excess material can easily be removed from both the buccal and palatal surfaces, so when the guide is removed, the mock-up is already clean and complete. For optimal flexibility and handling, the buccal area should maintain a thickness of around **4.5 mm**. With this approach, we can produce a printed guide that is both flexible and accurate, perfectly replicating what we used to achieve with a traditional silicone (PVS) index. The major difference lies in control: digitally, we can fine-tune the buccal thickness and make virtual cuts in the occlusal and palatal areas, resulting in easier clean-up and more efficient use of material.



TECHNIQUE	SILICONE INDEX	PRINTED INDEX
Cost	€/ \$ 4-6	€/ \$ 2-5 (a full tray)
Production time	~3 h + guides	~3.5 h (depending on the layer height)
Possible mistakes	• Material selection • Tray selection to place the silicone for thickness control • Distortions and bubbles • Incorrect cuts of the silicone index	• Inadequate design of the guides • Failures in post-processing protocol
Standardization of the result	• Thickness control is not standardized • The quantity of material is not precisely controllable	• The results will be standardized, since what is designed in CAD will be obtained
Precision & surface resolution	• The quality of the index is completely dependent on the skills of the operator • It requires the extra step of printing the model, then handmade silicone PVS index manufacturing	• The quality of the index depends on the layer height selected and the thoroughness of the post-processing protocol • No need to print models
Dependence	• Knowledge, experience, and skill are necessary to correctly manufacture guides	• Need a good design of the guides, and a Form 4B ecosystem (printer and post-processing) and IBT Flex Resin

Protocol 3: Printed Mock-up Shell With Tooth-Colored Resins




An alternative to create the mock-up that utilizes 3D printing is to print the shell directly with a tooth-colored resin. The development of biocompatible resins enables this technique and speeds up the process.

Tools needed to do this protocol:

DIGITAL CAPTURING AND DESIGN	DIGITAL DESIGN	DIGITAL PRODUCTION	CLINICAL TOOLS
Intraoral scanner Digital camera	CAD software for doing the smile design + to design the shells	Form 4B 3D printer Resin Tank Resin Mixer Flex Build Platform Form Wash (2nd Generation) Form Cure (2nd Generation) Premium Teeth Resin OR BEGO™ VarseoSmile® TriniQ® Resin Accessories: IPA (99%) Protective equipment	Flowable composite Light curing device Probe Explorer tool Tweezers Gloves Black marker

Workflow Steps

FIRST APPOINTMENT: DATA ACQUISITION



Intraoral scanning
Photography
Videos

PLAN AND DESIGN



Design the wax-up on the patient's scan, create the shells, and export the STLs for the three sections.

3D PRINT AND POST-PROCESS THE PRINTED SHELLS



Print the shell using Premium Teeth Resin or BEGO™ VarseoSmile® TriniQ® Resin.

DELIVERY



Use flowable composite to set the printed shell in place.

1. Design Considerations

In this case, we are focusing on producing the shell for the motivational mock-up. The shell is only an additive design. The process begins with a design split into three separate STL files, representing distinct sections of the arch (anterior, and two posteriors). Why three? Traditional full-arch shells (like in the early DSD workflows) often cracked or failed upon insertion, especially when arch anatomy varied. By dividing the shell into three sections, fit is far more reliable — no model, no index, just print and place.



Picture by Dr. Jesús Creagh

2. Material Selection

There are two tooth-colored resins available from Formlabs Dental to print using Form 4B: [Premium Teeth Resin](#) and [BEGO™ VarseoSmile® TriniQ® Resin](#). The differences of these two materials can be found in the following table.



RESIN	PREMIUM TEETH RESIN	BEGO™ VARSEOSMILE® TRINIQ® RESIN
Ideal for	Printing the shell for motivational or technical mock-up	
Other uses of the material	A nano-ceramic filled biocompatible material with enhanced aesthetics and translucency. Indicated for denture teeth , temporary full-arch implant-supported restorations (All-on-X appliances), temporary single units (crowns, inlays, onlays, and veneers), and and up to seven-unit temporary bridges.	A versatile ceramic-filled biocompatible material, indicated for temporary and permanent single units (crowns, inlays, onlays, and veneers) and bridges (up to three units for permanent indications), and denture teeth.
	Class II (US), Class IIa (EU)	Class II (US), Class IIa (EU).




RESIN	PREMIUM TEETH RESIN	BEGO™ VARSEOSMILE® TRINIQ® RESIN
Shades	Bleach A2, A3, B1	A2, A3, B1
Cost of cartridge	\$/€ 550/kg	\$/€ 599 / 750 g
Cost per part	\$/€ <1	\$/€ 1.20-2
Layer height recommended and estimated print time for three shells (one full arch)	Motivational mock-up: 100 µm, 21 min* Technical mock-up: 50 µm, 57 min* <i>*Please note that print times may vary depending on the geometry and height of the parts</i>	Motivational mock-up: 100 µm, 30 min* Technical mock-up: 50 µm, 1 h 33 min*
Post-processing time	Follow the Manufacturing Guide Wash: • Form Wash (2nd Generation): Wash for 10 min or until clean Cure: • Form Cure (2nd Generation): 5 min at 80 °C. OR • Form Cure (1st Generation)/Form Cure L: 30 min at 80 °C	Follow the official recommendations Wash: • Wash for 3 min and spray with fresh IPA Cure: • Form Cure (2nd Generation): 10 min at 60 °C OR • Form Cure (1st Generation)/Form Cure L: 2x 20 min at 60 °C
What type of mock-up can I use this for?	Motivational Mock-up: print at 100 µm Technical Mock-up: print at 50 µm	
Guidance for choosing the right resin	If you want to print mock-ups and indications, such as All-on-X appliances, provisional restorations, or teeth for dentures, then this resin would be a good fit. If you are looking for a material with high translucency, this would be the best fit. Dental professionals should verify regional availability of this material here .	If you are looking into printing permanent restorations as well as provisional ones, then this resin would be a good choice. Users must use the Flex Build Platform to print with this material. Dental professionals should verify regional availability of this material here .



For printing the shells, Premium Teeth Resin is the recommended choice as it provides translucency and natural aesthetics, especially in the B1 shade. The results out of the printer are already impressive, with a nice surface finish that polishes nicely with minimal effort.

3. Nesting in PreForm Dental

In this case, the shells need to be oriented so that their incisal edges and cusps face towards the build platform. Then, supports can be automatically generated using **Support > Auto-Generate All**.



Confirm that you don't have any support structures in undesired areas. You can remove these manually under the **Supports** menu by clicking **Edit** and selecting the support you want to remove.



4. 3D Printing and Post-Processing

Following nesting in PreForm Dental, the parts are printed and post-processed. The 3D printing workflow involves several steps: printer setup, printing, washing, drying, and post-curing. Total production time varies depending on the layer height chosen for printing.

4.1 Setting up Form 4B

Insert a build platform, resin tank, resin mixer, and Premium Teeth Resin cartridge in the printer. Once the printer is set up, send the print job from PreForm Dental to the printer.

When printing with biocompatible resins, a dedicated build platform, resin tank, and resin mixer are needed to ensure biocompatibility.

4.2 Post-processing

After the print job is completed, the next step to follow is post-processing which includes removing the parts from the build platform, washing, drying, and post-curing. Post-processing of parts printed with biocompatible resins requires following the steps as established in the [Premium Teeth Manufacturing Guide](#).


Wash:

- After parts are removed from the build platform, they need to be washed in IPA (99%) in the Form Wash (2nd Generation) at High Speed for 10 minutes or until clean. It's important to use a dedicated Form Wash for washing parts printed with Premium Teeth Resin. This Form Wash should not be shared with other materials and the exact timing presented in the [Manufacturing Guide](#) should be followed.
- Immediately after washing, dry the printed parts with compressed air to ensure that there is no remaining IPA or resin on the printed part, especially in the intaglio surface. This is critical to ensure proper fit of the shell during the clinical execution.
- Allow parts to air dry at room temperature for at least 30 minutes.


Cure:

- Parts need to be post-cured in the Form Cure (1st Generation) or Form Cure (2nd Generation) for the specific time and temperature as established in the [Manufacturing Guide](#). This final step will ensure parts are ready to be used clinically.

The following graphics summarize the overall process and final production time.



Technical Mock-up: Premium Teeth Resin at 50 μm



4.3 Part Finishing

4.3.1 Support removal:


Remove support structures with a cutting disc or very carefully use a pair of high-quality precision nippers to trim off support tips as close to the surface as possible. Grind down the touch points with a dual-cut carbide bur.





5. Clinical Execution

After polishing, the parts are ready to be used clinically.



No bonding or etching are needed. Just position the printed shells and use **flowable composite** to close the gaps.



Place the **anterior shell** and remove any excess material.



Light cure that section.



Next, position one of the **lateral shells**.



Remove excess material and light cure.



Finally, place the shell on the remaining side.



Remove excess material and light cure.



Once the shell is positioned and secured, continue with the photographic and video documentation.



Despite the **minimal thickness (1 mm)**, the **material strength, translucency, and low opacity** provide a highly convincing look that patients instantly connect with.

Once the documentation is completed with the motivational mock-up, the printed shell sections can be easily removed with the explorer tool.

This technique is fast to execute, precise, and the most efficient of all three protocols.

Conclusion

3D printing has removed the barriers that once made mock-up fabrication complex, time-consuming, or inconsistent. With the ability to choose from three streamlined protocols, dental professionals now have full control over materials, time, and workflow — whether it's selecting the right model resin, or deciding to skip bis-acryl altogether in favor of a high-precision 3D printed shell. The technology is here, it's accessible, and it's proven. Today, there's no reason not to leverage digital tools to deliver fast, clean, motivational and technical mock-ups that elevate patient experience and clinical outcomes. The choice is yours — customize your approach, and let 3D printing do the heavy lifting.



Get in Touch

Whether you need to produce occlusal splints in-house, or are looking for high-throughput production of dental models, we're here to help. Get in touch with a dental expert to find the right solution for your business.

[Find a Reseller](#)

Product Demo: Form 4B

Blazing speed meets unmatched accuracy in the next generation of dental 3D printing. Explore the latest innovations in our on-demand webinar.

[Watch the Webinar](#)

Request a Free Sample Part

See and feel Formlabs quality firsthand. We'll ship a free 3D printed sample part printed on the Form 4B to your office.

[Request a Free Sample Part](#)